



# Double-stranded RNA (dsRNA) ELISA kit (J2 based)

## Features

- Highly sensitive detection of dsRNA
- Detects plant and animal viruses
- Distinguishes bacterial from viral pathogens

## Ordering Information

### Catalog Number + (size)

10613002 (200 tests)

10613005 (500 tests)

### Format

ELISA Kit

### Species Reactivity

Ubiquitous

### Company Information

Exalpha Biologicals, Inc.

tel: 978.425.1370

info@exalpha.com

www.exalpha.com



## **Other Nucleic Acid Detection Kits & Reagents Available from Exalpha Biologicals**

Double-stranded RNA Detection Kit (K1 based)

10623002 (200 Tests)

10623005 (500 Tests)

Mouse anti double-stranded RNA (J2, J5 and K1) Set

10040200 (3x100 µg)

DNA Fragmentation Detection Kit

X2044K1 (30 Tests)

X2044K2 (60 Tests)

BrdU Cell Proliferation Assay Kit

X1327K1 (200 Tests)

X1327K2 (1000 Tests)

X1327K3 (5000 Tests)

BrdU Chemiluminescent Cell Proliferation Assay Kit

X1623K1 (200 Tests)

X1623K2 (1000 Tests)

X1623K3 (5000 Tests)

BrdU Immunohistochemistry Kit

X1545K.1 (50 Sections)



## Table of Contents

Intended Use	4
Storage of Kit Components	4
Background and Principle of the Assay	4-6
References	6
Materials Provided	6
Materials Required but Not Provided	6-7
Preparation of Reagents	7
Assay Protocol	7-10
Precautions and Recommendations	10
Ordering & Contact Information	11





### **Intended Use**

The Exalpha Biologicals, Inc. Double-stranded RNA (dsRNA) ELISA kit (J2 based) can be used to detect viral dsRNAs or large natural or synthetic dsRNAs of non-viral origin in nucleic acid extracts, as well as to detect the presence of undesired dsRNA molecules in artificially synthesized (m)RNA preparations. By using serial dilutions of the Poly (I:C) dsRNA standard (included in the kit) for calibration, quantitative estimates can be made.

**This assay is for research use only and not for use in diagnostic or therapeutic procedures.**

### **Storage of Kit Components**

Exalpha's Double-stranded RNA (dsRNA) ELISA Kit components are shipped on cold packs. Upon receipt, store entire kit at -20°C. Once the kit is thawed, you may keep it at 4°C for 5 days. For long-term storage, it is recommended to aliquot and freeze the antibodies and dsRNA component at -20°C.

### **Background and Principle of the Assay**

Based on the use of two double-stranded RNA (dsRNA)-specific monoclonal antibodies the dsRNA Detection Kit allows sensitive and selective detection of dsRNA molecules (larger than 30-40 bp), independent of their nucleotide composition and sequence. The detection is highly specific: dsRNA can be detected in nucleic acid extracts in the presence of 1.000-10.000-fold excess of other nucleic acids. This assay works on the sandwich-ELISA principle and uses the J2 (IgG2a) mouse monoclonal antibody to dsRNA as a catcher antibody. The monoclonal antibody K2 (IgM) is used as the detector antibody.



Over the past decade our double-stranded RNA (dsRNA) antibodies have been used extensively to detect and characterize plant and animal viruses with dsRNA genomes or intermediates. In addition, the anti-dsRNA antibodies can be used as a tool to detect pathogens, including detection in paraffin-embedded fixed tissue samples (Richardson et al. 2010).

The J2 anti-dsRNA IgG2a monoclonal antibody (Schönborn et al. 1991) has become the gold standard in dsRNA detection. It was used initially for the study of plant viruses, but since the seminal paper of Weber et al. in 2006, where J2 was used to show that all the positive strand RNA viruses tested produced copious amounts of dsRNA in infected cells, this antibody has been used extensively in a wide range of systems, as documented in over 200 scientific publications.

J2 can be used to detect dsRNA intermediates of viruses as diverse as Hepatitis C virus, Dengue virus, rhinovirus, Chikungunya virus, Rabies virus, Polio virus, Classic swine fever virus, Brome mosaic virus and many more in cultured cells and also in fixed paraffin-embedded histological samples.

J2 has been used to elucidate how anti-viral responses are initiated, what counter-strategies viruses have adopted to avoid them, and to explore the viral life cycle by enabling ultrastructural localization studies of viral nucleic acid replication sites (Knoops et al., 2011). J2 has also been recommended as a tool to detect whether an unknown pathogen is bacterial or viral in nature (Richardson et al., 2010). Recently J2 has also been used to monitor the removal of dsRNA from in vitro synthesised mRNA preparations that may have potential use in gene therapy



(Kariko et al., 2011). J2 has been used successfully in various immunocapture methods, such as ELISA.

## References

- J. Schönborn et al. Nucleic Acids Res.(1991)19, 2993.  
 S. J. Richardson et al. J Clin Virol. (2010) 49(3); 180.  
 F. Weber et al. J Virol (2006), 80(10):5059-64.  
 K. Knoops et al. J Virol. (2012) 86(5); 2474.  
 K. Karikó et al. Nucleic Acids Res. (2011) 39(21) e142.

## Materials Provided

Double-stranded RNA (dsRNA) ELISA kit (J2 based) is provided in 200 and 500 test size. Volumes listed below are for the 200 test kit followed by the 500 test kit.

	Component name and information	Part #	200 test kit	500 test kit
1.	Coating antibody	J0170	63 µl	156 µl
2.	Poly (I:C) dsRNA as positive control in RNase/DNase-free, sterile STE buffer. The concentration is 1µg/µl. (Store at -20°C or -80°C)	J0171	4 µl	10 µl
3.	dsRNA-specific detecting antibody (in RPMI + 5% FBS)	J0172	22 ml	55 ml
4.	HRP-conjugated F(ab') <sub>2</sub> Fragment of goat-anti mouse secondary antibody	J0173	3 µl	5 µl
5.	TMB substrate solution (store at +4°C, keep in dark)	J0005	22 ml	55 ml



## Materials Required But Not Provided

1. 2 ELISA plates (96 wells; e.g. Nunc Immunoplate F96 Maxisorp or Costar cat nr 2595)
2. Microtiter plate reader spectrophotometer with wavelength capability at 450 nm.
3. Single channel pipettes 10  $\mu$ l and 200  $\mu$ l.
4. Multichannel pipettes 200  $\mu$ l or squirt bottle.
5. Antigen (standard and sample) diluent (STE Buffer: 0.1 M NaCl, 1 mM EDTA, 50 mM Tris-HCl, pH 7.0)
6. Washing Buffer (PBS + 0.5% Tween 20; PBS: 10 mM Pi-buffer, pH 7.2, 0.15 M NaCl).
7. Secondary antibody dilution buffer (PBS+1% BSA).
8. Blocking buffer (PBS containing 1% BSA and 0.2% NaN<sub>3</sub>)
9. Optional storage buffer (PBS containing 0.2% NaN<sub>3</sub>)
10. Incubator allowing incubation at 37 °C.
11. 2 M H<sub>2</sub>SO<sub>4</sub>.

## Preparation of Reagents

1. Use DEPC treated MilliQ water to prepare STE (when applicable for your own sample preparation)
2. Sterilize PBS and STE by autoclaving or filter through a 0.2 micron filter
3. Prepare PBS+1% BSA, ELISA washing buffer and blocking buffer.
4. Prepare storage buffer if needed.

## Assay Protocol

1. Coating the plates
  - a. For two plates, add 63  $\mu$ l of J0170 into 21 ml PBS, mix well and immediately distribute 100  $\mu$ l/well in 2 ELISA plates. For five plates, add 156  $\mu$ l of J0170 into 52 ml PBS, mix well and immediately distribute 100  $\mu$ l/well in 5 ELISA plates.



- b. Cover the plates and incubate them overnight at 4 °C.
      - c. Discard contents of wells into waste. Add 100 µl/well 1% BSA in PBS + 0.2% NaN<sub>3</sub> to each well and incubate at 37 °C for 2 h to saturate any remaining free binding sites on the plate.
      - d. Discard the solution and wash plates 3 times with PBS + 0.5% Tween 20.
      - e. The plates can then be used directly or stored. For storage fill the wells with 200 µl/well PBS containing 0.2% sodium azide.
2. Preparation of the Poly (I:C) dsRNA Positive Control.
  - Prepare 1:3 serial dilutions from J0171 by using RNase/DNase-free, sterile STE buffer.
  - The dilution series of the dsRNA standard should be in the range of expected dsRNA concentration of your sample.
  - We propose starting with 30 ng dsRNA/well as the highest concentration and diluting down to below 0.01 ng dsRNA/well.
  - Dilutions should be freshly made for each assay.
3. Preparations of Sample
  - Prepare dilutions of your sample in STE (when necessary).
  - Cap and vortex all diluted standards and samples.
4. Wash Step
  - Remove the plastic foil from the ELISA plate
  - Discard contents of wells into waste. Wash plate 4 times with PBS + 0.5 % Tween 20 adding 250 µl washing solution/well. Discard the solution.
  - Do not allow wells to dry before adding the next solution.
5. Addition of the Antigen





- Transfer 100  $\mu$ l antigen or diluted standard to duplicated wells in the plate.
  - Cover and Incubate 1 hour at 37 °C.
6. Wash Step
- Discard contents of wells into waste. Wash plate 4 times with PBS + 0.5 % Tween 20 adding 250  $\mu$ l washing solution/well.
  - Do not allow wells to dry before adding the next solution.
7. Detector Antibody Addition
- Pipette 100  $\mu$ l undiluted J0172 into all wells.
  - Incubate 1 hour at 37 °C.
8. Dilute secondary antibody
- During the incubation (step 5) dilute J0173 by pipetting 1.3  $\mu$ l into 21 ml PBS + 1% BSA (no azide!) for 200 tests or 3.25  $\mu$ l J0173 into 52 ml PBS + 1% BSA (no azide!) for 500 tests.
9. Wash Step
- Discard contents of wells into waste. Wash plate 4 times with PBS + 0.5 % Tween 20 adding 250  $\mu$ l washing solution/well.
  - Do not allow wells to dry before adding the next solution
10. Secondary Antibody Addition
- Add 100  $\mu$ l diluted J0173 into each well.
  - Incubate 1 hour at 37 °C.
11. Final Wash Step
- Discard contents of wells into waste. Wash plate 4 times with PBS + 0.5 % Tween 20 adding 250  $\mu$ l washing solution/well.
  - Do not allow wells to dry before adding the next solution
  - Take care to remove all washing fluid after the last wash.
12. Development
- Add 100  $\mu$ l of J0005 into each well.  
Incubate for 5-60 minutes at room temperature in the dark.



### 13. Stop

- When the absorbance has reached the optimum level stop reaction by adding 100  $\mu$ l of 2M H<sub>2</sub>SO<sub>4</sub> to all wells.

### 14. Read

- Read absorbance at 450 nm, blanking on the zero standard.

## Precautions and Recommendations

1. All standards and samples should be assayed at least in duplicate.
2. Use clean, RNase-free micro-centrifuge tubes with cap.
3. Do not use buffers which contain NaN<sub>3</sub> as it will interfere with the final detection step.
4. Do not expose reagents to excessive light.
5. Wear disposable gloves and eye protection.
6. Do not use the kit beyond the expiration date.
7. Do not mix reagents from different kits.
8. Do not mouth pipette or ingest any of the reagents.
9. The buffers and reagents used in this kit contain anti-microbial and anti-fungal reagents. Care should be taken to prevent direct contact with these products.
10. Do not smoke, eat, or drink when performing the assay or in areas where samples or reagents are handled.
11. Human samples may be contaminated with infectious agents. Do not ingest, expose to open wounds, or breathe aerosols. Dispose of samples properly.
12. After completion of step 1: Cell Plating the plates can be stored without any loss of activity for one month. To store wrap plates in plastic foil and store them refrigerated at 4 °C. When stored plates are used, they must be thoroughly washed with PBS to remove all traces of NaN<sub>3</sub>.



## Ordering information

Catalog Number	Size
10613002	200 tests
10613005	500 tests

## Contact Information


Exalpha Biological, Inc.  
2 Shaker Road, Bldg. B101  
Shirley, MA 01464

*Telephone:* 978.425.1370

*Website:* [www.exalpha.com](http://www.exalpha.com)

*E-mail:* [info@exalpha.com](mailto:info@exalpha.com)

*Version 2.0*



Exalpha Biologicals,  
Inc. | 2 Shaker Rd.  
Bldg. B101 | Shirley,  
MA 01464  
Tel: 978.425.1370 |  
[info@exalpha.com](mailto:info@exalpha.com)  
[www.exalpha.com](http://www.exalpha.com)